

# USF CONNECT

## Shimadzu High Performance Liquid Chromatography (HPLC) LC20-AT



The HPLC is used to determine the amount of a chemical compound within a mixture of other chemicals, for example the exact amount of active pharmaceutical ingredient in a drug.

### Utilization

Install Date: November, 2007

Number of Samples Run: 900 (as of 4/1/08)

Frequency of Use: Weekly

A Data Station was purchased in conjunction with the HPLC.

This system includes a liquid chromatograph, degasser, autosampler, UV/Vis detector and CBM.